



*Royal
vibe*



QUICK START USER MANUAL

Table of Contents

ROYALVIBE KIT CONTENTS	3
QUICK START STEPS	4
Cellphone Requirements:	5
ROYALVIBE USER MANUAL	6
FIRST TIME USE / LOGIN	6
MAIN MENU	7
PAIRING YOUR DEVICE	8
INFORMATION MENU	10
TREATMENTS MENU	11
TREATMENT LIST SCREEN	12
TREATMENT LIST SCREEN continued	13
PRODUCT TERMS AND CONDITIONS:	14
INVESTMENT GUARANTEE:	14
DISCLAIMER:	14
MEDICAL INFORMATION & HISTORY FORM:	15
WARRANTEE:	16
CONTACT INFORMATION	17
LOG-IN DETAILS	17

ROYALVIBE KIT CONTENTS

- RoyalVibe
- Sound Health Decoder
- EnviroVibe
- Power Cable
- Connecting Cable (Gold Tip)
- Antenna (Red)
- 2 Therapy Balls
- 2 Connecting Cables – Blue (Negative) & Red(Positive)
- Encryption Key – Sound Health Decoder
- USB Cable – Sound Health Decoder

QUICK START STEPS

1. Plug the Blue RoyalVibe box into a power source via the power cable.



2. Connect the Antenna or EnviroVibe to the RoyalVibe with the Gold-tipped connecting cable.



3. Connect the therapy balls with the blue and red cables to the RoyalVibe (blue) box.



Red cable to the top (+/Positive) and Blue cable to the bottom (-/Negative)



4. Make sure the power switch on the RoyalVibe is switched on.



5. Volume knob adjusts the amplitude (strength). The higher the volume is set the stronger the signal sent would be.



NB: Therapy balls must not be in direct contact with each other when powered. This will cause a short-circuit

6. On your Android device make sure the Bluetooth is switched on and paired to device called RoyalVibe.
7. Download the RoyalVibe Application from the Google Play Store. It's **FREE** to download.
8. Register on www.HealthAdviceHub.com for your account.
 - You can register a maximum of 5 dependents under your profile, each having their own treatments and schedules.
9. On your RoyalVibe application, log into your Health Advice Hub account on the application. Find your treatment under the correct dependent and run the program.
10. If your treatments have been scheduled simply ensure that all devices are on, and programs will run. Scheduled treatments will over-ride any current programs running.
11. Preferably keep a dedicated Android Smartphone for RoyalVibe use only.

***NOTE:**


The therapy balls are for electro-therapy and are to be used with direct-contact on the affected area for treatment of muscles / joint relief. The red antenna is for all the sound frequency treatments

Cellphone Requirements:

- 1) Android 5 and higher
- 2) Bluetooth 4
- 3) 1 Gig Ram
- 4) 5 or 5 inch screen
- 5) Data Sim Card

ROYALVIBE USER MANUAL


FIRST TIME USE / LOGIN



- For first time registration, click on **REGISTER** button or login using **GUEST LOGIN**. The registration will take you to our website – **HEALTH ADVICE HUB**, where you will complete the process before coming back to **LOGIN**
- **The Health Advice Hub, CellQuicken and Terms & Conditions** buttons all take you to the relevant website links
- To log in click the **LOGIN** button




- Complete the user name and password (same as used on **HEALTH ADVICE HUB**) and click the **SUBMIT** button
- If you forgot your password, you can click on the **Forgot Password** button, which will take you to the website to reset your password.
- When login is successful, select the dependent and click **SELECT** button



MAIN MENU

(Swipe in from the left)




- The **TREATMENTS** button takes you to where you can choose between Public or Personal groupings or the Quick Treatment search.
- **SCHEDULED TREATMENTS** button takes you to your scheduled treatment list allowing you to rerun a treatment already done or that is yet to be done (it will only be treatments that are scheduled)
- **CHANGE SCHEDULE** button takes you to the website where you can change any of your already set scheduled treatments. You can also set other treatment groupings that you have already created.

- **INFORMATION** button takes you to the information menu
- **DEPENDENT:** Lets you choose / change the dependent
- **TUTORIAL:** Takes you various usage tutorials for download.
- **LOG USES:** Lets you see you log usage



➤ **SETTINGS:** Lets you change the flash rate, connect your Bluetooth, add your personalized message as indication that a program has ended and also choose language options.



➤ **SUPPORT HELP CHAT:** Allows for online Application support

➤ **SHARE:** Lets you share the application with others via different platforms.

➤ **VIDEOS:** Here you can view videos of interest and application assistance.

➤ **SIGNOUT:** Lets you sign out of the application to sign in as a different user.


PAIRING YOUR DEVICE

➤ Click on **BLUETOOTH** test to pair your app to your RoyalVibe box.

Your RoyalVibe needs to be switched on for this

➤ Click on the 1st **PAIR** button

You might be prompted to allow permissions for the app.



- After it scans for a while you will see a device RVxxxx

(Matching the serial number marked on the RoyalVibe box)




- Click on PAIR. The button will change as seen below. Click the back arrow to go back to the previous menu.



- Blue tick shows that the devices have been paired

- Now click on the 2nd Pair Button. Where you will see the below screen. Click on scan devices.




- Now pair the RVA with your device number.



INFORMATION MENU

- The **USER MANUAL** button will take you to the online manual. A PDF version is also available for download from the website.
- **HEALTH ADVICE HUB**, this will take you to the website where you can get advice and helpful information on a range of medical topics.
- **CELLQUICKEN** will take you to the CellQuicken website to see other products and services.
- **TERMS & CONDITIONS** takes you to the terms and conditions online.
- **SETTINGS** are for syncing your Bluetooth with the RoyalVibe Unit (RVU). This is also where you can change the channel on your device if you have multiple RVUs in your household.

This allows that 2 or more RVUs can be used simultaneously from different devices (one device per RVU at a time) in the same room without interfering with each other's unit. An RVU cannot be paired to multiple devices.




TREATMENTS MENU

When you have clicked on the **TREATMENTS** button on the main menu you will get to this page.

- ❖ **PUBLIC GROUPINGS** displays a list of treatments that is available for everyone to use
- ❖ **PERSONAL GROUPINGS** displays a list of groupings that you have created but not shared publicly
- ❖ **QUICK TREATMENT SEARCH** displays a list of all the available treatments allowing you to Add keywords and do a quick search.

POP UP ALERT:


Reminder that you must log on to an active internet connection and run a program live at least once every 14 days to ensure that you have the latest Data updates.



TREATMENT LIST SCREEN


When you have clicked on the **PUBLIC** or **PERSONAL GROUPINGS** button on the **SELECT TREATMENT LIST** screen you will land to this page. This will be a treatment list relevant to your choice in the previous screen and whether it is Treatments, Audio or Visual

- ❖ Here you can do a search or scroll down the list. To select one click on load
- ❖ You can add up to 5 treatments in the que to run consequetively, or just run one at a time.
- ❖ Once you have selected the treatment you want to run, you click on the **RUN** button at the bottom of the screen
- ❖ Here you can also choose to **RUN ONCE**, **LOOP** or **NO** (disabled) on your audio and visual treatments.
- ❖ When you are done with your selection, remember to press **RUN**.
- ❖ Once you have clicked the **RUN** button, you will see a pop up as follows:



TREATMENT LIST SCREEN continued

- ❖ This pop up screen as in the example on the right, will disappear and your program will start to run
- ❖ When you want to pause a running treatment, click the **PAUSE** button.
- ❖ To resume simply click the **RESUME** button and the countdown will continue.
- ❖ When pausing, your RoyalVibe light will stay lit. To stop a running treatment press **STOP**.
- ❖ The **PRE-RUN SELECTION** screen will appear, and the RoyalVibe light will switch off.



- ❖ At the top you'll see your Treatment name and running: 1 of 25 (for example)
- ❖ You will see a circle running with the remaining time left of your treatment displayed underneath, and, in the center, you'll see a larger timer showing the current running treatment's time left.
- ❖ Once it gets to 00:00 it will start running: 2 of 25 (for example) and the timer in the center will start counting down from a different time.



- ❖ The main counter at the bottom will continue to count down. Your **RoyalVibe** is currently running a treatment.

PRODUCT TERMS AND CONDITIONS:

CellQuicken manufactures the CellQuicken product range. While professional assistance is offered in using the devices (billable per service), the product range is designed as a tool for self-care.

As the same with conventional medicine, lifestyle adjustments are important and therefore, we encourage and endorse healthy living habits.

INVESTMENT GUARANTEE:

We know that cell functions adjust to new vibrational indoctrination over a period of time and therefore the user must use the device on a daily basis to have noticeable results. Some results are more evident within shorter periods than others. If a client does not feel any benefit after continuous use, they must contact the offices to ensure that they are using it correctly and that the treatment program is sufficient that the correct medical information has been provided.

DISCLAIMER:

Under no circumstances are we suggesting or trying to replace your medicine or therapy or trying to convince you that this alternate modality will cure you. No pregnant persons, persons who have had vital organ transplants or persons who cannot work with good knowledge, skills and insights to equipment and their health may operate or apply electro or sound therapy and/or use the device.

The EnviroVibe contains a magnetic field and any persons with a pacemaker or other implants need to be cautioned against use and remain a safe distance away.

CellQuicken cannot guarantee that the use of the CellQuicken product range is beneficial to a person's specific health condition or that you will be cured of any ailment.

By using the CellQuicken product range and making use of frequency generation systems you, the user, understands and accepts that you have no demand or expectation of curing any ailment. The client also understands that the use of the CellQuicken product range affords no negative physical or mental effects on the human body and that neither CellQuicken, nor their representatives will be held liable for such anomalies. You the user understands that that the Sound Health Decoder device is not a medical diagnosis but merely detects any imbalances within the range of testing points. Moreover, the user voluntarily accepts all responsibility for the use and application of all frequencies generated by our systems. The user understands that any product purchased is subject to CellQuicken's return policy as set out. All returns are subject to obtaining an authorization number from customer service via email to info@cellquicken.com

MEDICAL INFORMATION & HISTORY

FORM:

When requesting a treatment program; the information and programs written are not intended to be a substitute for professional medical advice, diagnosis or treatment. Never disregard professional medical advice, or delay in seeking it because of information or programs received from CellQuicken or any of its associates.

If you think you may be suffering from any medical condition, you should seek immediate medical attention. CellQuicken or any of their associate employees are not responsible or liable for any advice, course, treatment, diagnosis, program or any information that you obtain from us. Any verbal or written communication is to be seen as a comment and must be confirmed by your own physician.

You hereby agree that the frequency program written on your behalf is based on information given by you and that it is only to assist with treatment programs and in no way, constitutes a diagnosis or treatment protocol.

WARRANTEE:

The base unit of the CellQuicken RoyalVibe has a 5-year Warranty

The CellQuicken Analyser has a 1-year Warranty

This is a carry-in warranty and any cost for transport will be at the clients' account. CellQuicken also takes no responsibility for any damages incurred because of courier transporting and damages will be paid by the client.

Exclusions: Any damage due to power spike, dropping physical damage and or other neglect and manhandling.

VR glasses carry the original manufacturer warranty and/or return policies. CellQuicken and its subsidiaries are not liable for operations and/or maintenance of said equipment and including operation of Cellphone required to run product application.

CONTACT INFORMATION

Website: www.Cellquicken.com

Office Tel: 012 – 516 0095

Mobile: 062 634 0038

General Queries: service@Cellquicken.com

Medical Queries: client@Cellquicken.com

Technical Assistance: support@Cellquicken.com

Physical Address: Fintech Campus Bldg, 81 Botterklapper Street
Pretoria East, 0081, Pretoria, South Africa

Office Hours: Monday to Thursday 08:00 to 16:00
Friday 08:00 to 13:00

LOG-IN DETAILS

Website: HealthAdviceHub.com

Usage: RoyalVibe Profile

Username: _____

Password: _____

Gmail Account (Set up by Cellquicken):

Usage: RoyalVibe Application on Playstore
_____. RoyalVibe@Gmail.com

Password: _____